


8 种亚硝胺类化合物的分析

组份: N-亚硝基二甲胺、N-亚硝基甲基乙基胺、N-亚硝基二乙胺、N-亚硝基二正丙胺、N-亚硝基二正丁胺、N-亚硝基哌啶、N-亚硝基吗啉

色谱柱: MEGA-WAX 25m x 0.32mm 1.2um df

订货号: S-WAX-032-120-25



since
1980

MEGA s.n.c. - Capillary Columns Laboratory - Via Plinio, 29 - 20025 Legnano (MI) - Italy - Tel./Fax: (+39) 0331.54.79.24 - info@mega.mi.it

www.mega.mi.it

